# KIERO Plast

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Español

### Descripción

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## English

### Product description

**Kiero Plast** is a light-curing modelling composite paste for dental laboratory use. It is suitable for milling, interlocking and object transfers thanks to its pasty form. The modelling is cold made and therefore there is no thermal contraction. Cured objects can be tested in mouth. The blue colour guarantees a fast polymerization in depth and represents a good contrast with the model. Kiero Plast doesn't contract during the polymerization process due to its low monomer content. It burns out completely and therefore it can also be used in pressed ceramic. It cures with standard 320-500nm lamps and it can be combined with ordinary waxes once it's polymerized.

This ductile material can be applied universally; nevertheless it is suitable for pontics, bars and large implants structures and even when using pantographs. Kiero Plast means a perfect complement with the liquid KIEROForm composite when modelling in dental laboratory.

### Composition

Composition of acrylic resin and initiators.

### Areas of application

- Modelling of individual implant abutments and bars
- Modelling of splints for interlocking
- Modelling of telescopic crown
- Interlocking of objects
- Modelling of crowns and bridges for pantographic transmission

### Properties and advantages

**Kiero Plast** patterns are highly stable and precise and can be removed without distortion. Patterns made of light cured or standard waxes can be combined ideally. The material is immediately position ally stable, i.e. it can be applied quickly, sparingly and accurately. **Kiero Plast** burns out without residue even in thicker layers and guarantees precise casting results.

### Instructions for use

After removing the lid from the container take the needed amount with your fingers or an adequate instrument. After using close the lid again in order to avoid light curing of not used material. You can employ Kiero Plast easily because of it's putty consistence.

### The perfect surface

**Important:** Do not remove the inhibition layer from the surface; it provides an optimum bond between the layers and prevents stresses in the modelling resin. The inhibition layer should only be removed from the surface using alcohol or acetone after contouring of the pattern is complete or if wax is to be added to the pattern.

**To obtain optimum casting results, complete removing of the inhibition layer is essential. On large volumes patters object should be complemented with a thin wax coating.**

### Polymerization

Kiero plast polymerases with standard light machines between 320 and 500 nm. It polymerases without vacuum. Its polymerization times are:

- 60-90 sec. with stroboscopy light
- 2-3 min. with halogen light
- 3,5-5 min. with UVA lamps

**Please note that curing times depend on the emission power of the lamps used in light curing devices. Correct maintenance of these devices as well as regular replacement of lamps account for a constant curing. Please also note the instructions of the device manufacturer.**

### Warning

Contains phosphine oxide, Tricyclodecan Dimethanol Diacrylat. Avoid any eye- or skin contact with uncured material as well as with the inhibited layer. Wear suitable protective gloves, protective clothing and eye protection. Do not breath dust/ vapours. Switch on aspirator when grinding. Uncured Material is dangerous for the environment. Please pay attention to the valid material safety data sheet also.

### Contraindications, Interactions

The use of the product is contraindicated if known allergy to any ingredient exists. Interactions are not determined.

### Secondary effects

Product may cause allergic reaction.

### Storage Instructions

Store at 4-25°C.

Protect Kiero Plast container from heat and sunlight.

Do not use after the expiry date.

### Disposal instructions

Dispose of contents/ container in accordance with local regulation.

### Delivery unit

Kiero Plast: 30 g in light protected container

This product was developed for use in dentistry and should be used according to the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also personally responsible for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

### Legend

Lot number

Recommended storage temperature

Expiry date

Manufacturer

Read instructions for use

Exclamation mark, Warning

Protect from sunlight

Environment

### Fabricante / Hersteller / Manufacturer

Kuss Dental, S.L.  
c/ Isabel Colbrand, 10 Nave 95  
28050 Madrid - Spain  
Tel. + 34 91 738 23 17  
Fax + 34 91 738 23 18  
info@kuss-dental.com  
www.kuss-dental.com  
www.kuss-dental-shop.com

Vers.3 - MT-KP-V-17/02/2016