KIERO Block

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Español

Descripción

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English

Product description

Kiero Block is a one red component light-curing resin. It does not burn out. Kiero Block is easy to apply with accuracy on the model, because of its low viscosity. It is very economical because it can be applied straight from the cannula. It is well distinguishable because of its red colour. It can be cured in any light machine (320 – 500nm). It is hugely stable and also in high temperature processes, such as thermo-distortion. You can use any type of separator. The materials could only be applied on the model, but not in the mouth.

Composition

Formulation of acrylates and photoinitiators.

Applications

- It is ideal for modelling while using pantographs like ZirkonZahn, Ceramill, etc., for its good stability
- To cover pieces before sandblasting
- To fill in the pores in models or other corrections
- To close saw cuts
- All kind of blockings (dies, shoulders, etc)
- Spacer for producing soft flexible bleaching splints

Instructions for use

Kiero Block is worked over a dry and clean model. It can be used with any kind of separator. The material is applied straight from the syringe. Once the material is used, it is not necessary to take off the cannula. The polymerization of layers over 1mm might cause an unfinished polymerization, an excess of temperature's release or tensions. Therefore, a perfect adjustment of the piece cannot be guaranteed. To reduce dilatations while polymerization, it must be done with maximum 1 mm thick layers and cure them.

The inhibition layer:

Alcohol or acetone can be used to remove the inhibition layer or you can work out the modelation. Do not remove the inhibition layer from the surface until the work is done. It guarantees an optimum bond between the layers and avoids the creation of more model tensions.

Polymerization time:

Please keep in mind that polymerization timing depends precisely on the intensity of the light emission of the used lamp. The right maintenance of the machine, including the regular change of light bulbs, helps to get a good polymerization. Check the manufacturer's instructions.

Observations

Kiero Block contains acrylic resins which might create irritation and hyper-sensitivity. Avoid contact of uncured product or inhibition layer with eyes and mucus. It should not stay in the mouth for long. To avoid health problems, it is required to work only in suitable places equipped with vacuum cleaner. Please, follow rigorously the indications of the safety data sheets.

Contraindications

Sensibilities in people with intolerances to some of its components cannot be excluded. In those cases, stop using the product.

Storage

Store between 4 and 25 °C.
Protect the syringes of heat and direct sunlight.
Do not use after expiry date.

Pack contents

Box of 3 x 3 g syringes Kiero Block.

This product was developed for use in dentistry and should be used according to the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also personally responsible for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

Legend

	Lot number		Manufacturer
	Expiry date		Xi, irritant
	Read instructions for use		Recommended storage temperature

Fabricante/Hersteller/Manufacturer

Kuss Dental, S.L.
c/ Isabel Colbrand, 10 Nave 147
28050 Madrid - Spain
Tel. +34 91 736 23 17
Fax + 34 91 736 23 18
info@kuss-dental.com
www.kuss-dental.com

Vers.2 – 09/12/2014 – MT-KB-EAI