KIERO Form

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Español

Descripción

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English

Product description

Kiero Form is a light-curing one-component modelling composite gel for dental laboratory use. It is suitable for milling, interlocking and object transfers thanks to its pasty form. The modelling is cold made and therefore there is no thermal contraction. The blue colour guarantees a fast polymerization in depth and shows a good contrast with the model. It burns out completely and therefore it can also be used in pressed ceramic. **Kiero Form** cures with standard 320-500nm lamps and it can be combined with ordinary waxes once it is polymerized. The **Kiero Form** gel means a perfect complement with the **Kiero Plast** composite paste when modelling in dental laboratory.

Composition

Composition of acrylic resin and initiators.

Application area

- Modelling of attachments, telescopic or cone crowns and individual abutments
- Objects interlocking
- Modelling of adhesive bridges, Inlays and Onlays
- Modelling of crowns and bridges for pantographic transmission like ZirkonZahn, Ceramill, .
- All kind of blockings (dies, shoulders, retentive zones, etc)
- Spacer for producing soft flexible bleaching splints

Properties and advantages

Kiero Form patterns are highly stable and precise and can be removed without distortion. Patterns made of light cured or standard waxes can be combined and corrected ideally. The material is immediately positionally stable, i.e. it can be applied quickly, sparingly and accurately. **Kiero Form** burns out without residue even in thicker layers and guarantees precise casting results.

Instructions

Kiero Form is worked over a dry and clean model. It can be used with any kind of separator. The material is applied straight from the syringe. After taking the cap off, put the application tip on the syringe. It can only be used once. Close the syringe afterwards. The polymerization of layers over 1mm might cause an unfinished polymerization and excessive temperature's releases or tensions. Therefore, a perfect piece adjustment cannot be guaranteed. To reduce dilatations while polymerization, it must be done with maximum 1mm thick layers and cure them. Once the material is used, it is not necessary to take the cannula off.

Important: Do not remove the inhibition layer from the surface; it provides an optimum bond between the layers and prevents stresses in the modelling resin. The inhibition layer should only be removed from the surface by using alcohol or acetone after contouring the pattern is complete or if wax is to be added to the pattern. To obtain optimum casting results, complete removing of the inhibition layer is essential. On large volumes patters object should be complemented with a thin wax coating.

Polymerization

Kiero Form polymerizes with standard light machines between 320 and 500 nm. It polymerizes without vacuum. Its polymerization times are:

- 90 sec. with stroboscopy light
- 2 min. with halogen light
- 3,5-5 min. with UVA lamps

Please note that curing times depend on the used lamps emission power in light curing devices. Correct maintenance of these devices as well as regular replacement of lamps are necessary for a constant curing. Please also note the device manufacturer instructions.

Observations

Kiero Form contains acrylic resins which might create irritation and hypersensitivity. Avoid contact of uncured product or inhibition layer with eyes and mucus. It should not stay in the mouth for long. To avoid health problems, it is required to work only in suitable places equipped with vacuum cleaner. Please, follow carefully the indications of the safety data sheets. Harmful to aquatic life with long lasting effects.

Contraindications

Sensibilities in people with intolerances to some of its components cannot be excluded. In those cases, stop using the product.

Storage instructions

Store between 4 and 25 °C

Protect **Kiero Form** syringes of heat and direct sunlight

Do not use after expiry date.

Disposal note

Avoid release to the environment. Dispose of contents/container in accordance with local disposal requirements.

Delivery units

Box of 3 x 3 g. syringes
9 Cannulas

This product was developed for dentistry use and should be used according to the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also personally responsible for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

Legend

| | | | |
|---|---------------------------|---|---------------------------------|
|  | Lot number |  | Recommended storage temperature |
|  | Expiry date |  | Manufacturer |
|  | Read instructions for use |  | Protect against light |

Fabricante/Hersteller/Manufacturer

Kuss Dental, S.L.
c/ Isabel Colbrand, 10 Nave 147
28050 Madrid - Spain
Tel. +34 91 736 23 17
Fax + 34 91 736 23 18
info@kuss-dental.com
www.kuss-dental.com

Vers.2 - MT-KF-V-20150429