# DeltaSplint

Gebrauchsinformation

Deutsch

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Englisch

## Description

DELTA SPLINT is a MMA-free, light-curing modelling resin for the production of functional therapeutic occlusal splints. The ready-to-use DELTA SPLINT one-component material is characterized by its high accuracy of fit and dimensional stability. The uncomplicated and thus extremely fast fabrication of the splint is guaranteed by the composite-like consistency of the material. The stability of the material simplifies the modelling considerably and ensures the exact reproduction of the bite situation. The properties of DELTA Splint and its dosage form make it possible to produce occlusal splints in the shortest possible time.

## Composition

Oligomeric acrylates, inorganic and organic fillers, photoinitiators.

## Indication (Areas of Application):

- DELTA SPLINT is used for the production of functional therapeutic occlusal splints. It can be used independently or in combination with suitable thermoformed foils as a carrier material.
- It can also be used as a block-out material or for modelling any shape.

## Contra-indication

DELTA SPLINT must not be used if allergy to any of the ingredients is known.

## Side effects

If the irradiation with the polymerization device is too short or if the polymerization device is defective or unsuitable, the material has not hardened sufficiently and taste irritations may occur when wearing the splint.

## Application

### Processing without thermoforming splint

- Block out the plaster models with wax and sufficiently isolate them with DELTAALGIN alginate insulation.
- Lightly grease fingers, instruments and opposing bite with DELTALINA. This ensures comfortable and fast working.
- Apply the one-component material occlusally to the dental rim and model the splint in the articulator.
- Light cure the splint in the articulator for 10 minutes in the Metalight or another polymerization unit with UV-A tubes (350 - 400 nm).

- Remove the splint from the model and place it with the inside facing upwards. polymerize for another 10 minutes.
- Finish, grind and polish the splint as usual.
- As an alternative to polishing, the splint can be varnished with DELTAGLACE. Apply a thin layer of varnish with a brush and cure in Metalight for 10 minutes.

### Processing with thermoforming splint

- Roughen the prepared thermoforming foil by milling or sandblasting with 50µm or 110µ aluminium oxide - 2 bar pressure.
- Apply DELTABOND to the roughened surfaces with a brush.
- 2 minutes polymerize in Metalight (no inhibition layer)
- Model on DeltaSplint and light cure.

### Repair and rework of splints

- Both DeltaSplints plints and PMMA splints can be repaired.
- Roughen the surface and apply DELTABOND with a brush.
- Polymerise in Metalight for 2 minutes.
- Apply DELTA SPLINT and light cure.

## Warnings

DELTA SPLINT contains diphenyl(2,4,6-trimethyl benzoyl)phosphine oxide. May cause allergic reactions. Causes skin irritation. Harmful to aquatic life, with long-lasting effects. DELTABOND/DELTAGLACE contains methyl methacrylate 2-propenic acid, reaction products with penta erythrite, diphenyl (2,4,6-trimethyl benzoyl) phosphine oxide. May cause allergic skin reactions. Causes severe eye damage. Harmful to aquatic life, with long-lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. In order to prevent damage to health, only work at workplaces with extraction units, especially when milling and polishing.

## Processing information

For optimal processing, we recommend to grease the fingers lightly with Vaseline and to maintain a maximum processing temperature of 25°C. Powerful light devices with high polymerisation temperatures can negatively influence the accuracy of fit.

## Storage conditions

Store at 4-25°C.  
Protect DELTA SPLINT from heat and sunlight.  
Close the containers tightly after use (no sunlight).

## Shelf life

Do not use after expiration date.

## Delivery units

**Delta Splint** 12 bars à 9 g. One bar is sufficient for one splint.  
**DeltaBond** light-cured bonding material for repair and rework of DELTA SPLINT und

|                   |  |
|-------------------|--|
|                   | PMMA splints as well as for bond to thermoforming foils. (10 ml) |
| <b>DeltaGlaze</b> | light-cured varnish (10 ml)                                      |
| <b>DeltaAlgin</b> | plaster separating agent on alginate basis (50 g.)               |
| <b>DeltaLina</b>  | Separating agent for the hands (10 g.)                           |

The product was developed for use in dentistry and must be used according to the instructions for use. The Manufacturer does not accept any liability for damage caused by improper use of system and material. Furthermore, the user must independently check the material for its suitability and applicability for the intended purposes before using it, especially if he diverts from the process described in the instructions for use, or if his process is not described in the instructions for use.

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DeltaMed GmbH