ENGLISH

EPI-GUIDE, BIRESORBABLE BARRIER MATRIX
WITH OPLA® PATENTED ARCHITECTURE

INDICATIONS

Epi-Guide is indicated for use as an adjunct to periodontal restorative surgeries in the treatment of periodontal defects following established surgical procedures.

INTENDED USE:

Epi-Guide is a periodontal membrane that aids: regeneration of the periodontal attachment apparatus; augmentation around implants placed in immediate extraction sockets; augmentation of insufficient ridges for later implantation; and treatment of alveolar ridge deficiencies.

CONTRAINDICATIONS
Use of Epi-Guide is contraindicated in the presence of active infection where purulence is produced.

DESCRIPTION

Epi-Guide is a bioresorbable polymer synthesized from D,D-L,L-polylactic acid. Its unique, open-cell, patented architecture is a matrix of incomplete void spaces, randomly sized and shaped, which have multiple communications with neighboring void spaces.

The polymer of the Epi-Guide is of the metabolite lactic acid and is degraded to CO₂ and H₂O. The basic polymer has a long history of safe medical use. Epi-Guide is flexible, hydrophilic and initially exhibits a vigorous uptake of fluid blood when placed on the wound. Surface tension of the fluid blood that invests in the barrier dressing promotes and improves the adherence of the barrier/tooth surface interface.

HOW IT WORKS

Clinical investigations have revealed that, once implanted in the body, Epi-Guide provides support to the initial blood clot and maintain competence of collateral circulation. It provides a unique density-gradient three-dimensional construction designed to attract, trap and retain fibroblasts and epithelial cells while maintaining space around teeth for development of bone and periodontal support tissues.

Unlike passive or defensive barriers, Epi-Guide organizes fibroblasts and epithelial cells to help arrest epithelial migration during the later phases of healing. Epi-Guide maintains its architecture and structural integrity for 20 weeks after implantation with complete bioresorption occurring between 6-12 months.

OPENING THE PACKAGE

Open the package and gently remove the barrier using a sterile surgical tong. Place the device in the surgical sterile field in preparation for the procedure.

DIRECTIONS FOR USE

Preparation for Implant - After the patient has achieved anesthesia, a full-thickness mucoperiosteal flap is elevated. The flap should extend at least one tooth anterior and posterior to the defect(s). The defects should be thoroughly debrided of granulation tissue and the roots properly scaled and root-planed to remove calculus and diseased cementum.

After adequate site preparation, the Epi-Guide membrane material should be trimmed with scissors or a scalpel to adequately cover the defect. The material should extend approximately 2-3mm beyond the edge of the defect(s), mesially and distally as well as not normally.

Application: Once the material has been trimmed properly, orient the embossed side facing away from the tooth's surface and toward the gingiva, leaving the flat side of the membrane to closely approximate the tooth surface.

Stabilizing Suture Placement. Although employment of a stabilizing suture to secure the barrier is not necessary for most situations, the use and type of suture to is at the discretion of the clinician.

USE OF BONE GRAFT OR VOID FILLER MATERIAL

The use of bone grafts or bone graft substitutes in conjunction with Epi-Guide is at the discretion of the clinician.

SURGICAL POINTS TO REMEMBER

- thoroughly prepare surgical site
- trim membrane material to size
- orient embossed side facing toward the gingiva
- replace flaps
- provide patient with detailed instructions for post-operative care

POST-OPERATIVE REMINDERS

Patients are instructed to rinse with an anti-microbial mouth rinse between meals for two weeks post-surgery to control anaerobic bacteria. Patients should not brush or floss teeth in and around the treated area until the clinician feels sufficient healing has occurred.

Patients should return at one week post-surgery for evaluation and adjustments in oral hygiene instruction.

Patients should be monitored on a regular basis for evaluation, oral hygiene instructions, and coronal scaling and rubber cup polishing for the first eight weeks.

In the event of flap recession where the material becomes exposed, Epi-Guide will remain functional and should not be disturbed. In such cases, oral hygiene to control anaerobic bacteria should be intensified.

POSSIBLE COMPLICATIONS

Any periodontal surgery may result in the following complications: thermal sensitivity, gingival recession, flap sloughing, flap inflammation, resorption or alkylosis of the treated root, some loss of alveolar crest, perforation or abscess formation, infection, irregular gingival formation in the healing process. Complications associated with the use of anesthesia are possible.

PRODUCT CARE

Sterility: Epi-Guide is sterilized by exposure to gamma radiation, using a validated process to ensure a SAL of 10-6. **Do not resterilize.** Sterility is guaranteed unless blister package is opened or damaged.

Storage: Store Epi-Guide in a cool, dry location. Keep Epi-Guide away from heat, sunlight, office autoclave, and cabinets associated with hot lighting fixtures or hot plumbing. Do not use if temperature dot on package has changed from light gray or dark gray to black.

Caution: Single use only. Discard unused portions of the device.

Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury or illness.

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