KIEROLink

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Español

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English

Product description

KIEROLink is a light-curing bonding system for common dental alloys and titanium with denture acrylic. The system consists of Basic Paste, Opaquer powder pink, dark pink, A3, white and Opaquer Liquid.

Composition

Basic Paste: Methacrylate content with photoinitiator

Opaquer Powder: Inorganic filler with pigments and dibenzoyl-peroxide

Opaquer Liquid: Methylmethacrylate, acrylated oligomer, pentaerythritol tri-tetraacrylate

Applications

Improving the bond between the metal framework and denture acrylic when finishing removable CrCo frameworks and masking retainers with tooth-colored A3 opaquer.

Instructions for use

KIEROLink Basic Paste

Sandblast the metal surface to be covered with clean 50 µm or 110 µm corundum abrasive.

Then remove any dirt, e.g. dust or grease, with a steam cleaner or pure acetone.

Using a suitable brush apply one to two uniform layers of Basic Paste 100 µm – 300 µm (0.1 mm to 0.3 mm) thick to coat the dry, sandblasted and cleaned metal surface.

The paste is cured in a light-curing unit. High-output units with a stroboscope or halogen light cure the material in 3 minutes. Units with UVA tubes with a light range of 320 nm to 400 nm (e.g. Metalight, Kiero LPC, Kiero Light ECO) require 5 minutes.

KIEROLink Opaquer

Place Opaquer Liquid in a clean mixing vessel and mix with the powder to a viscous, spreadable consistency (recommended mixing ratio 1:1 by weight). Then apply the mixture uniformly to the metal surface coated with basic paste.

The opaquer is cured in a light-curing unit. The opaquer has a dual-curing system (light and heat). This ensures reliable curing even in undercut or shaded areas (stroboscope and halogen light: 5 minutes, UVA tubes: 10 minutes).

Preparing the metal surface in this way improves the bond between the metal framework and denture acrylic for finishing CrCo frameworks with removable restorations.

Note: Only use the mixture once. Do not mix dried opaquer again for reuse.

Risks and Safety instructions

KIEROLink should not be used if there are any signs of irritation or with a known allergy to methacrylate-based resins.

Do not inhale dust or methacrylate fumes produced during preparing or trimming. Avoid direct contact of the unpolymerized material with the skin, mucosa or eyes, as it could cause burns and sensitizations. Always use an extractor at the work area to prevent damage to health. Wear protective gloves and eye protection.

Causes damage to organs (lung) through prolonged or repeated exposure if inhaled. Do not breathe dust.

Storage

Store at 4-25°C.

Avoid direct sun light.

Always close containers after use and keep out of reach of children.

Do not use the material after the expiry date.

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This product was developed for use in dentistry and should be used according to the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also personally responsible for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

V1_26.05.21_MT_KL_E/I/A