# AcriFix

## Gebrauchsinformation

## Deutsch

### Produktbeschreibung

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### English

### Product description

*AcriFix* is a light-curing resin gel for use in dental offices and dental laboratories, particular on implant works, for the completion of small-sized saddles on partial dentures, for the fabrication of anterior immediate prostheses with a single missing tooth, for relining saddles, for repairing cracked or fractured acrylic dentures, for

fixing broken off teeth to full or partial dentures, for the quick repair of chipped anterior veneers on dentures, for repairing base parts or the correction of bubbles, pores and cracks on acrylic denture bases. *AcriFix* cures between 360 and 420 nm and can be cured using any light-curing unit working within quoted spectrum. *AcriFix* is cured without vacuum.

### Composition

Composition of acrylates and initiators.

### Areas of application

- Fabrication of anterior immediate prostheses with a single missing tooth
- Relining saddles of partial dentures
- Completion of small-sized saddles of partial or implant dentures
- Repairs caused by fractures or cracks on dentures or base parts
- Fixing of broken off teeth on full or partial dentures
- Quick repair of chipped anterior veneers on dentures
- Correction of bubbles, pores and cracks on acrylic denture bases

### Contra-indications

Patient and user: *AcriFix* should not be used in case of known allergies against one of the ingredients.

### Side effects

Patient and user: *AcriFix* may cause allergic reactions. User: Frequent and constant skin contact with unpolymerized Material may cause skin irritation.

### Interactions

Interactions are unknown.

### Properties and advantages

*AcriFix* enables to work quickly without using kneading silicon or having long-period curing times. The material is easy to apply and immediately positionally stable, i.e. it can be applied quickly, sparingly and accurately. *AcriFix* is free of MMA.

### Instructions for use

After removing the seal cap use the provided canulas for an accurate application. These are intended for single use only. After use seal the slide syringe again. Seek for a constant layering while applying the material (1-2 mm) . The curing time depends on the used light unit. Curing lights using stroboscope lamps as UniXS need 3 min, halogen lamp units as Spektramat 10 min and UVA-tube-based units (e.g. MetaLight) 15-20 min, depending on the number of tubes. Models made of white plaster support the polymerisation process. Polymerisation with hand-held lamps (halogen/plasma without filter) within the spectrum of 360 – 420 nm in the dental practice are possible as well.

**Notice:** Curing with LED light units is not possible.

**Attention:** When working on models, they shall be insulated with **AcriFix Separator**. When fixing acrylic parts with *AcriFix* as broken off teeth or fractured denture parts, **AcriFix Bonder** shall be used at any time before applying the resin. Apply the bonder with a brush on each acrylic part (or tooth) and light cure for approx. 3 min.

When fabricating immediate prostheses, teeth shall be sand-blasted and provided with chamfers to guarantee an adequate mechanical retention to the material. Afterwards **AcriFix Bonder** shall be applied with a brush on the teeth. Then cure as described above.

When relining a saddle, intermediate cure the material after applying *AcriFix*. Then remove the denture from the model and cure final as described above.

To apply material after curing, do not remove the inhibition layer on the surface. It guarantees an optimum bonding and avoids any stress to the resin. Upon completion remove the inhibition layer from the surface with rotary instruments or isopropanol.

### Warning

Pay attention to Safety Data Sheet.

*AcriFix* contains 1,4-Butandiol dimethacrylate, diphenyl(2,4,6-trimethylbenzoyl) phosphine oxide.

*AcriFix* Bonder contains methyl methacrylate 2-propenoic acid, reaction products with pentaerythrit, diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide.

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye damage. May cause respiratory irritation. Harmful to aquatic life with long lasting effects. Avoid breathing dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Gently wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention.

Use extraction system while working with rotary instruments on *AcriFix*.

#### Disposal

Dispose content/container in accordance with local disposal requirements.

#### Storage

Store at 4-25° C.  
Protect *AcriFix* syringes from heat and sunlight  
Do not use after the expiry date.

#### Pack contents

**AcriFix kit:**  
3 x 3 g *AcriFix* pink in light protected slide syringes  
5 ml *AcriFix* Bonder  
12 canulas

#### AcriFix accessories

3 g *AcriFix* A3 in light protected slide syringe  
3 g *AcriFix* Incisal in light protected slide syringe  
5 ml *AcriFix* Bonder  
50 ml *AcriFix* Separator

This product was developed for use in dentistry and should be used according to the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also personally responsible for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

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DeltaMed GmbH