

KIERO FORM

» Universal Light-Cure Modelling Composite «

Kiero Form is an universal light-cure composite gel for milling and attachment parts modelling and cast items blocking. Curing is carried out with high and low-power UVA light-emitting devices (320 y 500 nm). It can be applied directly from the syringe and easily lifted from the model. Once it's been light cured the composite becomes very hard and it can be worked out easily.

- ▶ **Kiero Form** modelling jobs are very stable and can be lifted up from de model without any deformation (test in mouth)
- ▶ Cured **Kiero Form** can be combined with commercially available waxes and Metacon
- ▶ **Kiero Form** burns out without leaving residues and it is also suitable for ceramic injection technology
- ▶ Very distinguishable because of its blue colour.

Indications

- ▶ Modelling of adhesive bridges, inlays and onlays
- ▶ Modelling of attachments, telescopic and conical crowns
- ▶ Intra- and extra-oral cast parts blocking (soldering)
- ▶ Metacon o **Kiero Plast** modelling corrections
- ▶ Defined gap elaboration for bleaching splints
- ▶ Retentive parts blocking
- ▶ Sensitive teeth covering while sandblasting
- ▶ Pores filling and plaster models saw cuts sealing



3 syringes of 3 gr.

For wider modellations we recommend **Kiero Plast**

