

ACRIFIX

» Always ready to use for acrylic prosthesis
in surgeries and laboratories «

AcriFix is an antiallergic light-curing gel for both dental surgeries and laboratories. It enables work to be performed quickly without prior walls or long polymerisation times and stands out due to its uncomplicated processing, high levels of stability and outstanding modelling properties. AcriFix does not contain any MMA, therefore it is particularly biocompatible. It can be cured easily in commercially available light-emitting devices (360-420 nm) without vacuum. Because of its compact form, it can be worked out easily and turns highly brilliant. It is available in pink, A3 and incisal colours.

Application areas

- ▶ Screw hole sealing in implant supported prosthesis
- ▶ Bubbles, pores, cracks or small defects corrections in denture bases
- ▶ Pink gums expansion in order to cover metal or little defects, specially in immediate load implants.
- ▶ Attachments and teeth fixing in mouth
- ▶ Bite corrections
- ▶ Teeth vertical dimension increase
- ▶ Production of partial dentures without prior walls
- ▶ Antiallergic prosthesis

Accessories

- ▶ 3 g Acrifix **pink** in a light-protected syringe
- ▶ 3 g Acrifix **A3** in a light-protected syringe
- ▶ 3 g Acrifix **incisor** in a light-protected syringe
- ▶ 5 ml bonder



3 syringes of 3 g with pink Acrifix
5 ml bonder
12 cannulas



Corrections



Implant supported
prosthesis



Partial prosthesis



Screw holes sealing in
pink or A3 colours



Protrusions



Exceptional
addings

Importer:



c/ Isabel Colbrand 10 Nave 95
28050 Madrid - Spain.
Tel. + 34 91 736 23 17
Fax. + 34 91 736 23 18
www.kuss-dental.com

Distributor:



Scan the QR to see the video